


Bacterial Genomic DNA Shearing

Example protocols and results are based on customer feedback.



Sonication Conditions:
Lane 1 – unsheared
Lane 2 – 5min Total Sonication
Lane 3 – 10min Total Sonication
*Samples run on 1% agarose gel, 100bp ladder.

Protocol Information

Cell Type: Bacterial Genomic DNA, 1ug/ml

Total Sample Volume: 500ul **Sonicator Amplitude Setting:** 60%

Sonication Pulse Rate: 10 seconds On, 10 seconds Off

Total Sonication On Time: 5min or 10min. Desired fragment size dependent. See resulting gel image above

Sample Process Temperature: 3°C

Customer notes:

- -Sonicated in 1.5ml mirocentrifuge tube polypropylene (commercially available).
- -Sample was in TE buffer
- -We use Illumina for all sequencing. NEB next gen kits work well with sonicated DNA using the Q800R.
- -System worked very well.