



ADJUSTABLE-VOLUME PIPETTE

pipetman[®] m

Progress made simple!



EASY TO USE

ADVANCED ERGONOMICS

TRUSTED RESULTS



GILSON[®]

EASY TO USE

ADVANCED ERGONOMICS

TRUSTED RESULTS



pipetman[®] m

Progress made simple!

PIPETMAN[®] M is the latest motorized member of the PIPETMAN family from Gilson. PIPETMAN M is a fusion between the renowned accuracy, precision and robustness of the PIPETMAN Classic and user-friendly functions that will simplify your work day. PIPETMAN M is so easy to use that you won't have to change the way you pipette, whether using the single or the multichannel version!

Ready to use

Start using it right out of the box!

1 Set the volume by rotating the push-button left or right.

2 Pipette with a simple press of the push-button.



3 Access additional features such as repetitive mode, mix function or easy readjustment with the mode-button.



Easy to use

No need to change how you pipette!

USE PIPETMAN M LIKE A MECHANICAL PIPETTE

Keep your technique consistent by using only one button to adjust the volume and to pipette.

ALL YOU NEED IN ONE PIPETTE TO OBTAIN TRUSTWORTHY RESULTS!



Simplify your daily routine with the intuitive software: no instruction necessary!

Easy navigation with direct access to all supplementary functions thanks to the mode button:

- 4 essential pipetting modes: pipetting, repetitive, mix and reverse.
- Piston speed adjustment.
- Configuration: Service, Easy readjustment, Maintenance.


Advanced ergonomics

Minimal effort, maximum comfort!

Discover the differences! PIPETMAN M offers the same legendary robustness, accuracy and precision as the PIPETMAN Classic with additional features that will make your work easier...

OPTIMIZED BALANCE


Lightweight, but built with strong materials and components make PIPETMAN M one of the lightest motorized pipettes ever. Its optimized weight gives you a perfectly balanced pipetting sensation.



IMPROVED COMFORT

PIPETMAN M contributes to lowering Repetitive Strain Injuries (RSI):

- Thanks to its motor, PIPETMAN M offers you pipetting with minimized pipetting forces.
- Low ejection and volume setting forces assist the progress of your daily work.



Always operational

Save time, improve efficiency!

USE WHILE CHARGING


PIPETMAN M can be used while charging without compromising performance.

LONGER USE BETWEEN CHARGE CYCLES

The Lithium-ion battery inside PIPETMAN M has been designed to ensure 900 pipetting cycles at maximum speed (up to seven 96 well plates) between charge cycles. PIPETMAN M charges 80 % of its full battery capacity in less than an hour and takes only three hours to fully charge.

Clever GLP features

Good Laboratory Practice (GLP) features ensure the high quality and reliability of your work.



pipetman[®]

More than a multichannel pipette... it's a PIPETMAN M!


Enjoy the ease of use of PIPETMAN M combined with the well-known advantages of the PIPETMAN Neo multichannel lower part!

1 Innovative patented piston drive mechanism

Absolute consistency in volume distributed!

The specific rack rail mechanism guarantees identical performance on each channel.

Thanks to an exclusive patented piston drive mechanism, the **pressure is applied consistently and simultaneously on each piston**, resolving the common downfall of multichannel pipettes: disparities between channels of the volume distributed.



2 Ingenious tip ejector technology

Best protection from Repetitive Strain Injuries (RSI)!

For more comfort, the ejector spacer has been designed to apply equal ejection forces over each tip sequentially.




This symmetric system **drastically reduces tip ejection forces and risks of RSI.**

3 Unique tip-holder design

No need to purchase a dedicated tip brand!

Successive rings have been positioned depending on the pipette volume range to guaranty the snug fit of tips throughout the protocols.



Complete range

Available in single and multichannel version!

8 models cover a range of volumes from 0.5 µL to 1000 µL.

Model	Ordering reference	Volume range	Color code
SINGLE CHANNEL MODELS			
P10M	F81001	0.5 µL to 10 µL	
P20M	F81002	2 µL to 20 µL	
P200M	F81003	20 µL to 200 µL	
P1000M	F81004	100 µL to 1000 µL	
MULTICHANNEL MODELS			
P8x20M	F81005	1 µL to 20 µL	
P12x20M	F81006	1 µL to 20 µL	
P8x200M	F81007	20 µL to 200 µL	
P12x200M	F81008	20 µL to 200 µL	



power carrousel™

PIPETMAN M exclusive charging stand!

This rotating stand, which holds up to 5 PIPETMAN M (single and multichannel), stores and charges your entire pipette set at one time on a minimum bench space!

Protect your PIPETMAN M from shocks and contaminations, and improve your efficiency by keeping your tools always operational!

Ordering reference: FB1001

Trusted results

Best results guaranteed with a performance exceeding ISO 8655 recommendations.

PIPETMAN M single and multi's specifications are more stringent than those required to be ISO 8655 compliant.

PIPETMAN M helps to avoid variable results between users, thanks to motorized pipetting.

This special feature guarantees highly reproducible results, especially on small volumes.

Model	Volume		Maximum Permissible Errors			
			Gilson		ISO 8655	
			Systematic error (µL)	Random error (µL)	Systematic error (µL)	Random error (µL)
SINGLE CHANNEL MODELS						
P10M	Min.:	0.5	± 0.040	≤ 0.013	± 0.12	≤ 0.08
		1	± 0.025	≤ 0.012	± 0.12	≤ 0.08
		5	± 0.060	≤ 0.020	± 0.12	≤ 0.08
	Max.:	10	± 0.080	≤ 0.025	± 0.12	≤ 0.08
P20M	Min.:	2	± 0.075	≤ 0.025	± 0.2	≤ 0.1
		10	± 0.100	≤ 0.035	± 0.2	≤ 0.1
	Max.:	20	± 0.150	≤ 0.050	± 0.2	≤ 0.1
P200M	Min.:	20	± 0.40	≤ 0.15	± 1.6	≤ 0.6
		100	± 0.80	≤ 0.22	± 1.6	≤ 0.6
	Max.:	200	± 1.00	≤ 0.26	± 1.6	≤ 0.6
P1000M	Min.:	100	± 2.5	≤ 0.4	± 8.0	≤ 3.0
		500	± 3.5	≤ 0.7	± 8.0	≤ 3.0
	Max.:	1000	± 5.0	≤ 1.1	± 8.0	≤ 3.0
MULTICHANNEL MODELS						
P8x20M	Min.:	1	± 0.08	≤ 0.05	± 0.4	≤ 0.2
		10	± 0.15	≤ 0.10	± 0.4	≤ 0.2
P12x20M	Max.:	20	± 0.25	≤ 0.12	± 0.4	≤ 0.2
P8x200M	Min.:	20	± 0.50	≤ 0.16	± 3.2	≤ 1.2
		100	± 1.00	≤ 0.30	± 3.2	≤ 1.2
	Max.:	200	± 2.00	≤ 0.50	± 3.2	≤ 1.2

Also available to you:

Utilizing the PIPETMAN M Stand Adapter is one more way to simplify your work day!

PIPETMAN M Stand Adapter has been designed to fit perfectly on Gilson CARROUSEL stand and SINGLE™ pipette holder (Ref. F161406). It combines two features in one: it holds and it charges your pipette at the same time!



Stand Adapter
(Ref. F807023)



CARROUSEL™ stand
(Ref. F161401)

Gilson's Pipetting System

Your protocols and experiments require more than just generic pipette tips. Indeed, tips are essential components of your pipetting system. According to ISO 8655-2*:

"The maximum permissible errors always apply to the total system of piston pipette and tip. When using pipette tips not delivered by the pipette supplier, the supplier declaration or certificate of conformity does not apply".

*ISO 8655 is an international standard for "piston operated apparatus". It defines, among other things, terms and definitions regarding piston operated pipettes and their testing, testing procedure and maximum permissible errors.

Pipette model	Recommended PIPETMAN DIAMOND tips		
	Model (scale 1/2)	Volume range	
SINGLE CHANNEL MODELS			
P10M	 D10	0.1 – 10 µL	
	 DL10	0.1 – 20 µL	
	 DF10	0.1 – 10 µL	
	 DFL10	0.1 – 10 µL	
P20M	 D200	2 – 200 µL	
	 DF30	2 – 30 µL	
P200M	 D200	2 – 200 µL	
	 D300	20 – 300 µL	
	 DF200	20 – 200 µL	
P1000M	 DF300	20 – 300 µL	
	 D1000	100 – 1000 µL	
	 D1200	100 – 1200 µL	
	 DF1000	100 – 1000 µL	
	 DF1200	100 – 1200 µL	
	 DF1200	100 – 1200 µL	
	MULTICHANNEL MODELS		
	P8x20M & P12x20M	 DL10	0.1 – 20 µL
 DFL10		0.1 – 10 µL	
 DF30		2 – 30 µL	
 D200		2 – 200 µL	
P8x200M & P12x200M	 D200	2 – 200 µL	
	 D300	20 – 300 µL	
	 DF100	10 – 100 µL	
	 DF200	20 – 200 µL	
	 DF300	20 – 300 µL	

As a member of our PIPETMAN® family,

 is compatible with most leading tip brands.

000 «Диаэм»

Москва
 ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
 +7 (812) 372-6040
 spb@dia-m.ru

Новосибирск
 +7 (383) 328-0048
 nsk@dia-m.ru

Воронеж
 +7 (473) 232-4412
 vrn@dia-m.ru

Йошкар-Ола
 +7 (927) 880-3676
 nba@dia-m.ru

Красноярск
 +7 (923) 303-0152
 krsk@dia-m.ru

Казань
 +7 (843) 210-2080
 kazan@dia-m.ru

Ростов-на-Дону
 +7 (863) 303-5500
 rnd@dia-m.ru

Екатеринбург
 +7 (912) 658-7606
 ekb@dia-m.ru

Кемерово
 +7 (923) 158-6753
 kemerovo@dia-m.ru

Армения
 +7 (094) 01-0173
 armenia@dia-m.ru

