


### Product Specification

Item Number	FS-E30
Product Name	PCR Film Sealer (Portable)
Dimensions(L*W*H)	47.2x11.9x14.3cm
Size/ctn	54*18*22cm
Weight	7kg
Compatibility	96-well/384-well PCR reaction plate, ELISA plate.
Product Features	<ul style="list-style-type: none"> <li>• Simple operation and reliable sealing</li> <li>• Automatic rolling mode, no need for personnel to operate</li> <li>• One time complete sealing within 30 seconds, high efficiency</li> <li>• Intelligent touch of entrance and exit control; Preventing problems such as samples not being placed properly or personnel making mistakes in operation, timely touch the switch can withdraw the samples from the working chamber</li> <li>• Small footprint, light weight, easy placing and moving at will</li> <li>• The status light can visually judge the machine status.</li> </ul> <p>Yellow light : standby state,  Green light - running state,  Red light - abnormal state,  Whenever red light is on, the machine can be restarted directly. And it can be used normally after automatic initialization.</p>
Storage Condition	Room Temperature
Brand	Servicebio
Image	

**000 «Диаэм»**
**Москва**

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ [sales@dia-m.ru](mailto:sales@dia-m.ru)
**[www.dia-m.ru](http://www.dia-m.ru)**
**С.-Петербург**  
+7 (812) 372-6040  
[spb@dia-m.ru](mailto:spb@dia-m.ru)
**Новосибирск**  
+7(383) 328-0048  
[nsk@dia-m.ru](mailto:nsk@dia-m.ru)
**Воронеж**  
+7 (473) 232-4412  
[vrn@dia-m.ru](mailto:vrn@dia-m.ru)
**Йошкар-Ола**  
+7 (927) 880-3676  
[nba@dia-m.ru](mailto:nba@dia-m.ru)
**Красноярск**  
+7(923) 303-0152  
[krsk@dia-m.ru](mailto:krsk@dia-m.ru)
**Казань**  
+7(843) 210-2080  
[kazan@dia-m.ru](mailto:kazan@dia-m.ru)
**Ростов-на-Дону**  
+7 (863) 303-5500  
[rnd@dia-m.ru](mailto:rnd@dia-m.ru)
**Екатеринбург**  
+7 (912) 658-7606  
[ekb@dia-m.ru](mailto:ekb@dia-m.ru)
**Кемерово**  
+7 (923) 158-6753  
[kemerovo@dia-m.ru](mailto:kemerovo@dia-m.ru)
**Армения**  
+7 (094) 01-0173  
[armenia@dia-m.ru](mailto:armenia@dia-m.ru)
