

Media	Description	Cat #	Type	Fill Vol	Packaging		
Microbial Enumeration Tests							
Soybean-Casein Digest Broth	Trypticase™ Soy Broth	299107	Bottle, serum	100mL	10 each		
		299416	Bottle, screw cap/septum	100mL	10 each		
		299113	Bottle, screw cap/septum	500mL	10 each		
		211768	Dehydrated	500g	1 each		
		211771	Dehydrated	5lb	1 each		
		211772	Dehydrated	25lb	1 each		
	Tryptic Soy Broth	257107	30mL Vial	20mL	50 each		
		257488	Bottle, wide mouth	90mL	10 each		
		257486	Bottle, wide mouth	100mL	10 each		
		211824	Dehydrated	100g	1 each		
		211825	Dehydrated	500g	1 each		
		211822	Dehydrated	2kg	1 each		
		211823	Dehydrated	10kg	1 each		
		Soybean-Casein Digest Agar	Trypticase™ Soy Agar	221185	Plate	100mm style	20 each
221283	Plate			100mm style	100 each		
222205	Plate, BD™ Sterile Pack			100mm style	100 each		
299099	Bottle, screw cap/septum			500mL	10 each		
211043	Dehydrated			500g	1 each		
211046	Dehydrated			5lb	1 each		
Tryptic Soy Agar	211047		Dehydrated	25lb	1 each		
	236940		Dehydrated	100g	1 each		
	236950		Dehydrated	500g	1 each		
	236920		Dehydrated	2kg	1 each		
	236930		Dehydrated	10kg	1 each		
	Potato Dextrose Agar		Potato Dextrose Agar	299906	Bottle, screw cap/septum	500mL	10 each
				213300	Dehydrated	100g	1 each
				213400	Dehydrated	500g	1 each
Buffers and Solutions							
Phosphate Buffer pH 7.2	Phosphate Buffer pH 7.2 (Stock Solution)	257385	Bottle, crimp/stopper	100mL	25 each		
	Phosphate Buffer pH 7.2	214973	Bottle, screw cap/septum	500mL	10 each		
Buffered Sodium Chloride-Peptone Solution pH 7.0	Buffered Sodium Chloride-Peptone Solution pH 7.0	257086	Bottle, crimp/stopper	100mL	25 each		
		257087	Bottle, crimp/stopper	500mL	10 each		
Isopropyl Myristate	Isopropyl Myristate	257098	Bottle, crimp/stopper	500mL	10 each		



Media Solutions for USP <61> and <62> Compliance

Microbiological Examination of Nonsterile Products

Microbiology – It's what we do.

Find out what we can do for you.



For more information on BD media solutions or to place an order, please contact your local BD sales manager or visit us on the web at www.bd.com/ds.

1 United States Pharmacopeial Convention, Inc. 2009. The United States Pharmacopeia 32/The National Formulary 27, 2009. The United States Pharmacopeial Convention, Rockville, MD.
 2 European Pharmacopeia, 6th Edition. European Directorate for the Quality of Medicine, Council of Europe, 226 Avenue de Colmar BP907, F-67029 Strasbourg Cedex 1, France.
 3 Japanese Pharmacopoeia, Fifteenth Edition.



BD Diagnostics
 7 Loveton Circle
 Sparks, MD 21152-0999
 800.638.8663
www.bd.com/ds

BD Dehydrated and Prepared Media for Nonsterile Product Testing



Helping all people live healthy lives

USP/EP/JP Compliance

In 2006, the United States Pharmacopeial Convention published a revised Chapter 61, and introduced a new Chapter 62, that covers the Microbiological Examination of Nonsterile Products. These USP chapters became effective in May 2009, and are harmonized with the European Pharmacopeia and the Japanese Pharmacopeia.

BD has reviewed and revised, where applicable, the labeling, formulations and performance testing for the media that are specified in these new chapters. We stand ready to assist you in complying with this new testing requirement as well as with other USP¹/EP²/JP³ harmonized requirements in the future.




Quality Assurance

Mixing culture media ingredients is an exacting and delicate task, similar to cooking a delicious gourmet meal. BD starts the process with rigid quality control over raw materials and ingredients. Next, BD uses specially designed machinery to manufacture and control the output of consistently high quality products. The experience of dedicated line operators, laboratory and research technologists, and over 180 years of combined Microbiology experience goes into each container of BD Media. Finally, our ISO 9000 quality systems and cGMP manufacturing practices provide the consistency and control that allows BD to deliver the highest quality media available. We are proud of our facilities and our products, and welcome our customers to audit our facilities.

New BD Wide Mouth Bottles now available!

- 100mL and 90mL TSB Wide Mouth Bottles with a validated Sterility Assurance Level (SAL) of 10⁻⁶
- A twist-off cap and an opening of 48mm allows for easier access



USP CHAPTERS <61> AND <62> BD MEDIA FORMULATIONS

Media	Description	Cat #	Type	Fill Vol	Packaging
Test for <i>Escherichia coli</i>					
MacConkey Broth	MacConkey Broth	215177	Bottle	100mL	10 each
		220100	Dehydrated	500g	1 each
MacConkey Agar	MacConkey Agar	215197	Plate	100mm style	20 each
		297064	Plate	100mm style	100 each
		212123	Dehydrated	500g	1 each
		212122	Dehydrated	2kg	1 each
		275300	Dehydrated	10kg	1 each
		211387	Dehydrated	500g	1 each
		211390	Dehydrated	5lb	1 each
		211391	Dehydrated	25lb	1 each
Test for Bile-tolerant Gram-negative Bacteria					
Violet Red Bile Glucose Agar	Violet Red Bile Glucose Agar	215198	Plate	100mm style	20 each
		215053	Plate	100mm style	100 each
		218661	Dehydrated	500g	1 each
Mossel Enterobacteriaceae Enrichment Broth	EE Broth Mossel	292627	Bottle, wide mouth	90mL	10 each
		256620	Dehydrated	500g	1 each
		297005	Dehydrated	500g	1 each
Test for <i>Salmonella</i> spp					
Rappaport Vassiliadis Salmonella Enrichment Broth	RVS Soy Broth	215199	K Tube	10mL	10 each
		214943	Dehydrated	500g	1 each
Xylose-Lysine-Deoxycholate Agar	XLD Agar	221192	Plate	100mm style	20 each
		221284	Plate	100mm style	100 each
		278850	Dehydrated	500g	1 each
		278820	Dehydrated	2kg	1 each
		278830	Dehydrated	10kg	1 each
Test for <i>Pseudomonas aeruginosa</i>					
Cetrimide Agar	Pseudosel™ Agar Cetrimide Agar Base	297882	Plate	100mm style	10 each
		285420	Dehydrated	500g	1 each
Test for <i>Staphylococcus aureus</i>					
Mannitol Salt Agar	Mannitol Salt Agar	221173	Plate	100mm style	20 each
		221271	Plate	100mm style	100 each
		211407	Dehydrated	500g	1 each
		211410	Dehydrated	5lb	1 each
		293689	Dehydrated	25lb	1 each
Test for Clostridia					
Reinforced Medium for Clostridia	Reinforced Medium for Clostridia	215192	Bottle, screw cap/septum	100mL	10 each
		218081	Dehydrated	500g	1 each
Columbia Agar	Columbia Agar	215191	Plate	100mm style	20 each
		211124	Dehydrated	500g	1 each
		211125	Dehydrated	5lb	1 each
		211126	Dehydrated	25lb	1 each
Test for <i>Candida albicans</i>					
Sabouraud Dextrose Broth	Sabouraud Dextrose Broth	215193	Bottle, screw cap/septum	100mL	10 each
		238220	Dehydrated	100g	1 each
		238230	Dehydrated	500g	1 each
		238210	Dehydrated	2kg	1 each
Sabouraud Dextrose Agar	Sabouraud Dextrose Agar	221180	Plate	100mm style	20 each
		221278	Plate	100mm style	100 each
		292336	Plate, BD™ Sterile Pack	100mm style	100 each
		211584	Dehydrated	500g	1 each
		211585	Dehydrated	5lb	1 each
		210940	Dehydrated	100g	1 each
		210950	Dehydrated	500g	1 each
		211661	Dehydrated	2kg	1 each
		210930	Dehydrated	10kg	1 each

continued on reverse