

Compatible Chart of Filter Tips and Pipettes

Pipette Brand	Pipette Channel	Pipette Specification	Pipette SN/ Model	Specification of Servicebio Filter Tips						
				10μL Short (filter)	10μL Ext. (filter)	20μL (filter)	100μL (filter)	200μL (filter)	1000μL Short (filter)	1000μL Ext. (filter)
Thermo Finpipette	Single-Channel	1-10μL	4640000	√	√					
		2-20μL	4640020		√					
		10-100μL	4640040			√	√			
		20-200μL	4640050			√	√	√		
		100-1000μL	4640060						√	√
	Multi- Channel	1-10μL	4662000	√	√					
		10-100μL	4662020			√	√			
	30-300μL	4662030				√	√			
Brand	Single-Channel	0.1-2.5μL	705869 D-2.5		√					
		0.5-10μL	705870 D-10	√	√					
		2-20μL	705872 D-20			√				
		5-50μL	705873 D-50			√				
		10-100μL	705874 D-100			√				
		20-200μL	705878 D-200			√	√	√		
	100-1000μL	705880 D-1000						√	√	
	Multi- Channel	0.5-10μL	705900 M8-10	√	√					
		10-100μL	705908 M8-100			√	√			
		30-300μL	705912 M8-300				√	√		
Eppendorf	Single-Channel	0.1-2.5μL	3120000216		√					
		0.5-10μL	3120000224	√	√					
		2-20μL	3120000232			√				
		10-100μL	3120000240			√	√			
		20-200μL	3120000259			√	√	√		
		100-1000μL	3120000267						√	√
	Multi- Channel	0.5-10μL	3122000019		√					
		10-100μL	3122000035			√	√			
		30-300μL	3122000051					√		
Gilson Pipeetman	Single-Channel	0.2-2μL	P2	√	√					
		1-10μL	P10	√	√					
		2-20μL	P20			√				
		10-100μL	P100			√				
		20-200μL	P200			√	√	√		
		100-1000μL	P1000						√	√
	Multi- Channel	0.5-10μL	P10	√	√					
		2-20μL	P20			√				
	20-200μL	P200			√	√	√			
Rainin	Single-Channel	0.5-10μL	SL-10 PL+	√	√					
		2-20μL	SL-20 PL+			√				
		10-100μL	SL-100 PL+			√				
		20-200μL	SL-200 PL+			√	√	√		
		100-1000μL	SL-1000 PL+						√	√
	Single-Channel- XLS	0.5-10μL	SL-10 XLS	√	√					
		20-200μL	SL-200 XLS			√	√	√		
	100-1000μL	SL-1000 XLSC						√	√	
DragonLab	Single-Channel	0.5-10μL	7010101004	√	√					
		2-20μL	7010101005			√				
		5-50ul	7010101006			√				
		10-100μL	7010101008			√	√			
		20-200μL	7010101009			√	√	√		
		100-1000μL	7010101014						√	√
	Multi- Channel	0.5-10μL	7010103004	√	√					
		5-50μL	7010103006			√				
	50-300μL	7010103012				√	√			
Socorex	Single-Channel	0.5-10μL	825.001	√	√					
		20-200μL	825.02			√	√	√		
		100-1000μL	825.1						√	√

Disclaimer:

- Servicebio pipette tips are compatible with pipettes brands and models been tested as the above list.
 - Servicebio does not warrant the suitability for any particular application pipettes which haven't been listed on the compatibility chart.
- Latest Update in Jan. 2022

Compatible Chart of Non-Filter Tips and Pipettes

Pipette Brand	Pipette Channel	Pipette Specification	Pipette SN/	Specification of Servicebio Non-Filter Tips						
			Model	10μL Short	10μL Ext.	20μL	100μL	200μL	1000μL Short	1000μL Ext.
Thermo Finpipette	Single Channel	1-10μL	4640000	√	√					
		2-20μL	4640020		√					
		10-100μL	4640040			√	√			
		20-200μL	4640050			√	√	√		
		100-1000μL	4640060						√	√
	Multi- Channel	1-10μL	4662000	√	√					
		10-100μL	4662020			√	√			
		30-300μL	4662030				√	√		
Brand	Single Channel	0.1-2.5μL	705869 D-2.5	√	√					
		0.5-10μL	705870 D-10	√	√					
		2-20μL	705872 D-20			√				
		5-50μL	705873 D-50			√				
		10-100μL	705874 D-100			√	√			
		20-200μL	705878 D-200			√	√	√		
		100-1000μL	705880 D-1000						√	√
	Multi- Channel	0.5-10μL	705900 M8-10	√	√					
		10-100μL	705908 M8-100			√	√			
		30-300μL	705912 M8-300				√	√		
Eppendorf	Single Channel	0.1-2.5μL	3120000216	√	√					
		0.5-10μL	3120000224	√	√					
		2-20μL	3120000232			√				
		10-100μL	3120000240			√	√			
		20-200μL	3120000259			√	√	√		
		100-1000μL	3120000267						√	√
	Multi- Channel	0.5-10μL	3122000019		√					
		10-100μL	3122000035			√	√			
		30-300μL	3122000051				√	√		
Gilson Pipetman	Single Channel	0.2-2μL	P2	√	√					
		1-10μL	P10	√	√					
		2-20μL	P20			√				
		10-100μL	P100			√	√			
		20-200μL	P200			√	√	√		
		100-1000μL	P1000						√	√
	Multi- Channel	0.5-10μL	P10	√	√					
		2-20μL	P20			√				
		20-200μL	P200			√	√	√		
Rainin	Single Channel	0.5-10μL	SL-10 PL+	√	√					
		2-20μL	SL-20 PL+			√				
		10-100μL	SL-100 PL+			√	√			
		20-200μL	SL-200 PL+			√	√	√		
		100-1000μL	SL-1000 PL+						√	√
	Single Channel- XLS	0.5-10μL	SL-10 XLS	√	√					
		20-200μL	SL-200 XLS			√	√	√		
		100-1000μL	SL-1000 XLSC						√	√
	DragonLab	Single Channel	0.5-10μL	7010101004	√	√				
2-20μL			7010101005			√				
5-50ul			7010101006			√				
10-100μL			7010101008			√	√			
20-200μL			7010101009			√	√	√		
100-1000μL			7010101014						√	√
Multi- Channel		0.5-10μL	7010103004	√	√					
		5-50μL	7010103006			√				
		50-300μL	7010103012				√	√		
Socorex	Single Channel	0.5-10μL	825.001	√	√					
		20-200μL	825.02			√	√	√		
		100-1000μL	825.1						√	√
Disclaimer:										
1.Servicebio pipette tips are compatible with pipettes brands and models been tested as the above list.										
2. Servicebio does not warrant the suitability for any particular application pipettes which haven't been listed on the compatibility chart.										
Latest Update in Jan, 2022										

Disclaimer:

- Servicebio pipette tips are compatible with pipettes brands and models been tested as the above list.
- Servicebio does not warrant the suitability for any particular application pipettes which haven't been listed on the compatibility chart.

Latest Update in Jan. 2022

000 «Диаэм»

 Москва
 ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru
С.-Петербург
 +7 (812) 372-6040
 spb@dia-m.ru

Новосибирск
 +7 (383) 328-0048
 nsk@dia-m.ru

Воронеж
 +7 (473) 232-4412
 vrn@dia-m.ru

Йошкар-Ола
 +7 (927) 880-3676
 nba@dia-m.ru

Красноярск
 +7 (923) 303-0152
 krsk@dia-m.ru

Казань
 +7 (843) 210-2080
 kazan@dia-m.ru

Ростов-на-Дону
 +7 (863) 303-5500
 rnd@dia-m.ru

Екатеринбург
 +7 (912) 658-7606
 ekb@dia-m.ru

Кемерово
 +7 (923) 158-6753
 kemerovo@dia-m.ru

Армения
 +7 (094) 01-0173
 armenia@dia-m.ru
