

Corning® Bar Coded Storage Tubes

A system solution for Automated Sample Management

ДИА-М
современная лаборатория

www.dia-m.ru
заказ on-line


Corning Life Sciences is proud to provide you with the Corning 1D/2D Bar Coded Storage tube line. Compared to the competition, the 1D/2D storage line is a superior storage solution designed to provide maximum identification.

Features include:

- ▶ **MAXIMUM IDENTIFICATION** – Synchronized 2D and linear bar code – along with marking spot
- ▶ **MAXIMUM INFORMATION** – 14 x 14 dot 2D bar code
- ▶ **BAR CODE STABILITY** – Laser-etched, not an attached label
- ▶ **COMPOUND COMPATIBILITY** – Polypropylene construction throughout, making it very inert

Product Information

	1.3 mL 1D/2D Tube	384 Round Tube	384 Square Tube
Tube material	Polypropylene	Polypropylene	Polypropylene
Temperature range	-180° to +110°C (screw cap) -80° to +110°C (TPE cap)	-80° to +110°C (cap)	-80° to +110°C (cap)
Tube length	43.9 mm	14.0 mm	14.0 mm
Tube height	52.7 mm with screw cap	16.9 mm with cap	-
	47.5 mm with TPE cap	17.2 mm total height with tubes, caps, no lid	17.5 mm total height with tubes, heat seal, no lid
	53.7 mm total height with tubes and caps, no lid	23.1 mm total height with tubes, caps and lid	23.1 mm total height with tubes, heat seal and lid
	58.3 mm total height with tubes, caps and lid		
Volume	1.3 mL	120 µL	130 µL
Working volume	1.0 mL	100 µL	110 µL
1D codes	Linear, human readable	–	–
2D code	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits	14 x 14 dots, 14 digits



CORNING

Corning® 96 1D/2D Bar Coded Storage Tubes Ordering Information

Cat. No.	Product Description	Qty/Pk	Pk/Cs	Qty/Cs
8500	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, without screw caps, bulk pack	96/pk	10 pk/cs	960 tubes/cs
8501	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, without screw caps, racked with cover	96/rack	10 racks/cs	960 tubes/cs
8502	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, with screw caps, bulk pack	96/pk	10 pk/cs	960 tubes/cs
8503	96 1D/2D Bar Coded Storage Tubes, 1.3 mL, with screw caps, racked with cover	96/rack	10 racks/cs	960 tubes/cs
8504	96 Storage Tube Screw Cap, polypropylene with O-Ring, bulk pack	96/pk	50 pk/cs	4800 caps/cs
8505	96 Storage Tube Screw Cap, polypropylene with O-Ring, on mat, 96 caps per mat, clear with cover	480/pk	10 pk/cs	4800 caps/cs
8506	Thermoplastic elastomer (TPE) 96 Storage Tube Septum Cap, bulk pack	96/pk	100 pk/cs	9600 caps/cs
8507	Thermoplastic elastomer (TPE) 96 Storage Tube Septum Cap on mat, 96 caps per mat	960/pk	5 pk/cs	4800 caps/cs
8508	Empty racks with lids for 96 Screw Capped or Thermoplastic Elastomers (TPE) Capped Storage Tubes	5 racks/pk	2 pk/cs	10 racks/cs
8509	384 2D Bar Coded Storage Tubes, round, without plug caps, racked	384/pk	20 pk/cs	7680 tubes/cs
8510	384 2D Bar Coded Storage Tubes, square, racked	384/pk	20 pk/cs	7680 tubes/cs
8511	384 Storage Tube Plug Cap for round tube, bulk pack	384/pk	100 pk/cs	38400 caps/cs
8512	384 Storage Tube Plug Cap on mats for round tube, 384 caps per mat	3840/pk	10 pk/cs	38400 caps/cs
8513	Empty racks with lids for 384 round or square Storage Tubes	1 rack/pk	20 pk/cs	20 racks/cs

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

