

Flask, Boiling, Round Bottom, Narrow Neck

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Complies with ISO 1773
- With beaded rim

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
233.202.01	50	51	26	95	10
233.202.02	100	64	26	110	10
233.203.02*	150	74	28	118	10
233.202.03	250	85	34	143	10
233.202.04	500	105	34	168	10
233.202.05	1000	131	42	200	10
233.202.06	2000	166	50	240	6
233.202.07*	3000	185	50	260	1
233.202.08*	5000	223	50	305	1
233.202.09	10000	279	65	380	1
233.202.10*	20000	350	75	510	1

* Non DIN



Buy Now!



Flask, Boiling, Round Bottom, Wide Neck

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Complies with ISO 24450
- With beaded rim

Item Code	Capacity (ml)	Body Dia (mm)	Neck Dia (mm)	Height (mm)	Pack of
234.202.01	50	51	34	95	10
234.202.02	100	64	34	110	10
234.202.03	250	85	50	143	10
234.202.04	500	105	50	168	10
234.202.05	1000	131	50	200	10
234.202.07	2000	166	76	240	6



Buy Now!



Flasks



Buy Now!

Flask, Boiling, Round Bottom, Single Neck with Joint

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Complies with DIN / ISO 4797 & USP standard

Item Code	Capacity (ml)	Socket Size	Approx. Overall Height (mm)	Dia (mm)	Pack of
057.202.05	25	14/23	85	41	10
057.202.06*	25	19/26	85	41	10
057.202.07*	25	24/29	85	41	10
057.202.08	50	14/23	90	51	10
057.202.09	50	19/26	90	51	10
057.202.10	50	24/29	90	51	10
057.202.10A	50	29/32	90	51	10
057.202.11	100	14/23	105	64	10
057.202.12	100	19/26	105	64	10
057.202.13	100	24/29	105	64	10
057.202.14	100	29/32	105	64	10
057.202.16*	150	19/26	137	75	10
057.202.17*	150	24/29	137	75	10
057.202.18*	150	29/32	137	75	10
057.202.20*	250	14/23	140	85	10
057.202.21	250	19/26	140	85	10
057.202.22	250	24/29	140	85	10
057.202.23	250	29/32	140	85	10
057.202.24*	250	34/35	140	85	10
057.202.25*	500	19/26	163	105	10
057.202.26	500	24/29	163	105	10
057.202.27	500	29/32	163	105	10
057.202.28*	500	34/35	163	105	10
057.202.29*	1000	19/26	200	131	10
057.202.30	1000	24/29	200	131	10
057.202.31	1000	29/32	200	131	10
057.202.32*	1000	34/35	200	131	10
057.202.33	2000	24/29	240	165	6
057.202.34	2000	29/32	240	165	6
057.202.35*	2000	34/35	240	165	6
057.202.36A*	3000	24/29	260	185	1
057.202.36*	3000	29/32	260	185	1
057.202.37*	3000	45/40	260	185	1
057.202.40A*	5000	24/29	305	223	1
057.202.40	5000	29/32	305	223	1
057.202.40B	5000	34/35	305	223	1
057.202.41	5000	45/40	305	223	1
057.202.44*	10000	29/32	380	279	1
057.202.45A*	10000	34/35	380	279	1
057.202.45	10000	45/40	380	279	1
057.202.45B*	10000	55/44	380	279	1
057.202.47B*	20000	55/44	435	345	1

* Non DIN

Flask, Round Bottom, Single Neck with Joint

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- These Flasks come with Ground Joint manufactured as per ASTM E-1403 type 1 specifications
- The Flasks are printed in white colour enamel


Item Code	Capacity (ml)	Socket Size	Approx. Overall Height (mm)	Dia (mm)	Pack of
057.202.57	50	24/40	100	51	2
057.202.58	100	24/40	118	64	2
057.202.58A	200	24/40	128	75	2
057.202.59	250	24/40	140	85	2
057.202.59A	300	24/40	148	88	2
057.202.60	500	24/40	158	105	2
057.202.61	1000	24/40	185	130	1
057.202.62	2000	24/40	215	165	1
057.470.55B	25	24/40	90	41	2
057.470.09	50	19/22	80	51	2
057.470.11	100	14/20	100	64	2
057.470.12	100	19/22	100	64	2
057.470.20	250	14/20	120	85	2
057.470.21	250	19/22	120	85	2
057.470.25	500	19/22	140	105	2
057.470.25A	500	14/20	140	105	2
057.470.53	10	14/20	59	31	2
057.470.53A	10	19/22	59	31	2
057.470.55	25	14/20	75	41	2
057.470.55A	25	19/22	75	41	2
057.470.52	5	14/20	46	25	2



Flask, Round Bottom, 2 Neck Angular

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Complies with DIN 12394 & USP standard
- Available with two neck angular

Item Code	Capacity (ml)	Socket Size	Side Socket	Approx. Overall Height (mm)	Pack of
059.202.01*	50	24/29	14/23	92	1
059.202.02*	50	29/32	14/23	92	1
059.202.03*	100	24/29	14/23	107	1
059.202.04*	100	24/29	19/26	107	1
059.202.05*	100	29/32	14/23	107	1
059.202.06*	250	24/29	14/23	142	1
059.202.07*	250	24/29	19/26	142	1
059.202.08	250	29/32	14/23	142	1
059.202.09*	500	24/29	14/23	165	1
059.202.10*	500	24/29	19/26	165	1
059.202.11	500	29/32	14/23	165	1
059.202.12*	1000	24/29	14/23	203	1
059.202.13*	1000	24/29	19/26	203	1
059.202.14	1000	29/32	14/23	203	1
059.202.15A*	2000	24/29	19/26	243	1
059.202.15*	2000	34/35	19/26	243	1
059.202.20*	3000	34/35	19/26	260	1
059.202.20A*	3000	34/35	24/29	260	1
059.202.24*	5000	34/35	19/26	305	1
059.202.24A*	5000	34/35	24/29	305	1
059.202.28*	10000	34/35	24/29	380	1
059.202.29*	10000	55/44	24/29	380	1
059.202.30*	20000	55/44	24/29	435	1



PS : These flasks can also be provided in other joint sizes.

* Non DIN

Flasks

Flask, Round Bottom, 2 Necks Angular *New*

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Round bot tom flask with ASTM joint sizes
- Angular



Buy Now!



Item Code	Capacity (ml)	Socket Size	Side Socket	Pack of
059.470.00	50	14/20	14/20	1
059.470.00A	25	14/20	14/20	1
059.470.03	100	14/20	14/20	1
059.470.05A	100	24/40	24/40	1
059.470.08	250	24/40	24/40	1
059.470.10	500	24/40	24/40	1
059.470.15	1000	24/40	24/40	1

Flask, Round Bottom, 2 Neck Parallel

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Complies with DIN 12394
- Has Two Neck Parallel



Buy Now!



Item Code	Capacity (ml)	Socket Size	Side Socket	Approx. Overall Height (mm)	Pack of
059.212.00	50	24/29	14/23	90	1
059.212.00A	50	29/32	14/23	90	1
059.212.01*	100	24/29	14/23	105	1
059.212.01A*	100	24/29	19/26	105	1
059.212.01B*	100	29/32	14/23	105	1
059.212.02*	250	24/29	14/23	140	1
059.212.02A*	250	24/29	19/26	140	1
059.212.02B	250	29/32	14/23	140	1
059.212.03*	500	24/29	14/23	163	1
059.212.03A*	500	24/29	19/26	163	1
059.212.03B	500	29/32	14/23	163	1
059.212.04*	1000	24/29	14/23	200	1
059.212.04A*	1000	24/29	19/26	200	1
059.212.04B	1000	29/32	14/23	200	1
059.212.05*	2000	24/29	19/26	240	1
059.212.05A*	2000	34/35	19/26	240	1
059.212.06*	3000	34/35	19/26	260	1
059.212.06A*	3000	34/35	24/29	260	1
059.212.07*	5000	34/35	19/26	305	1
059.212.07A*	5000	34/35	24/29	305	1
059.212.08*	10000	34/35	24/29	380	1
059.212.08A*	10000	55/44	24/29	380	1
059.212.09*	20000	55/44	24/29	435	1

* Non DIN

Flask, Round Bottom, 3 Necks Angular

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Available with three neck angular
- Complies with DIN 12394

Item Code	Capacity (ml)	Centre Size	Side Socket	Approx. Overall Height (mm)	Pack of
060.202.01*	100	19/26	14/23	105	1
060.202.02*	100	24/29	14/23	105	1
060.202.03*	100	24/29	19/26	105	1
060.202.02A*	100	14/23	14/23	105	1
060.202.03A*	100	29/32	14/23	105	1
060.202.04*	250	19/26	19/26	140	1
060.202.04A*	250	19/26	14/23	140	1
060.202.05*	250	24/29	14/23	140	1
060.202.06*	250	24/29	19/26	140	1
060.202.07	250	29/32	14/23	140	1
060.202.08	250	29/32	19/26	140	1
060.202.08A*	250	29/32	29/32	140	1
060.202.09*	500	24/29	14/23	163	1
060.202.10*	500	24/29	19/26	163	1
060.202.11	500	29/32	14/23	163	1
060.202.12	500	29/32	19/26	163	1
060.202.12A*	500	29/32	29/32	163	1
060.202.13*	1000	24/29	14/23	200	1
060.202.14*	1000	24/29	19/26	200	1
060.202.15	1000	29/32	14/23	200	1
060.202.16	1000	29/32	19/26	200	1
060.202.17*	1000	29/32	29/32	200	1
060.202.18*	2000	24/29	19/26	240	1
060.202.19*	2000	29/32	24/29	240	1
060.202.20*	2000	29/32	29/32	240	1
060.202.26*	3000	34/35	19/26	260	1
060.202.26A*	3000	34/35	24/29	260	1
060.202.29B*	5000	34/35	19/26	305	1
060.202.29*	5000	34/35	24/29	305	1
060.202.32*	10000	34/35	24/29	380	1
060.202.33*	10000	45/40	24/29	380	1
060.202.34*	20000	55/44	24/29	435	1

* Non DIN



ASTM Sizes

Item Code	Capacity (ml)	Socket Size	Side Socket	Pack of
060.470.00	25	14/20	14/20	1
060.470.00A	50	14/20	14/20	1
060.470.01	100	14/20	14/20	1
060.470.01A	100	19/22	19/22	1
060.470.03	100	24/40	19/22	1
060.470.03A	100	24/40	24/40	1
060.470.04	250	14/20	14/20	1
060.470.04A	250	19/22	19/22	1
060.470.08A	250	24/40	24/40	1
060.470.09	500	24/40	24/40	1
060.470.10	500	19/22	19/22	1
060.470.19	1000	24/40	24/40	1



Flasks

Flask, Round Bottom, 3 Neck Parallel

- Made from ASTM E-438 Type 1, Boro 3.3 heat resistance glass
- Available Three Neck Parallel



Buy Now!

Item Code	Capacity (ml)	Centre Size	Side Socket	Approx. Overall Height (mm)	Pack of
061.202.00A*	100	19/26	14/23	105	1
061.202.00*	100	14/23	14/23	105	1
061.202.00B*	100	24/29	14/23	105	1
061.202.00C*	100	24/29	19/26	105	1
061.202.01*	250	24/29	19/26	140	1
061.202.01A*	250	19/26	19/26	140	1
061.202.02	250	29/32	14/23	140	1
061.202.02A*	250	29/32	29/32	140	1
061.202.03*	500	24/29	19/26	163	1
061.202.03A*	500	19/26	19/26	163	1
061.202.04	500	29/32	14/23	163	1
061.202.05	500	29/32	19/26	163	1
061.202.05A*	500	29/32	29/32	163	1
061.202.06*	1000	24/29	19/26	200	1
061.202.07	1000	29/32	14/23	200	1
061.202.08	1000	29/32	19/26	200	1
061.202.09*	1000	29/32	24/29	200	1
061.202.10*	1000	29/32	29/32	200	1
061.202.11*	2000	24/29	19/26	240	1
061.202.12	2000	29/32	19/26	240	1
061.202.13*	2000	29/32	29/32	240	1
061.202.14*	2000	34/35	19/26	240	1
061.202.15*	2000	34/35	24/29	240	1
061.202.17*	3000	34/35	19/26	260	1
061.202.18*	3000	34/35	24/29	260	1
061.202.19A*	5000	34/35	19/26	305	1
061.202.19*	5000	34/35	24/29	305	1
061.202.26*	10000	34/35	24/29	380	1
061.202.26A*	10000	45/40	24/29	380	1
061.202.27*	20000	55/44	24/29	435	1

* Non DIN

PS : All 2 Neck and 3 Neck Flasks can be provided in ASTM joint sizes on request

ASTM Sizes

Item Code	Capacity (ml)	Socket Size	Side Socket	Pack of
061.470.00	100	14/20	14/20	1
061.470.00A	100	24/40	24/40	1
061.402.06	250	24/40	24/40	1
061.402.09	500	24/40	24/40	1
061.470.05A	500	29/42	29/42	1
061.402.12	1000	24/40	24/40	1
061.470.10	1000	29/42	29/42	1
061.402.15	2000	24/40	24/40	1



Buy Now!

These flasks can also be provided in other joint sizes

Flask, Round Bottom, Four Necks Angular

- Complies with DIN 12392
- Centre Neck & Three Angled Side Necks
- With Interchangeable Joints
- Made From Heat Resistant, Low Expansion Borosilicate 3.3 Glass.

Item Code	Capacity (ml)	Interchangeable		Pack of
		Centre Neck	Side Neck	
060.240.00*	50	24 / 29	24 / 29	1
060.240.01*	100	19 / 26	14 / 23	1
060.240.02*	100	24 / 29	14 / 23	1
060.240.02A*	100	14 / 23	14 / 23	1
060.240.03*	100	24 / 29	19 / 26	1
060.240.03A*	100	29 / 32	14 / 23	1
060.240.04*	250	19 / 26	19 / 26	1
060.240.04A*	250	19 / 26	14 / 23	1
060.240.05*	250	24 / 29	14 / 23	1
060.240.06*	250	24 / 29	19 / 26	1
060.240.07*	250	29 / 32	14 / 23	1
060.240.08*	250	29 / 32	19 / 26	1
060.240.08A*	250	29 / 32	29 / 32	1
060.240.09*	500	24 / 29	14 / 23	1
060.240.10*	500	24 / 29	19 / 26	1
060.240.10A*	500	24 / 29	24 / 29	1
060.240.11*	500	29 / 32	14 / 23	1
060.240.12	500	29 / 32	19 / 26	1
060.240.12A	500	29 / 32	29 / 32	1
060.240.12B*	500	34 / 35	24 / 29	1
060.240.13*	1000	24 / 29	14 / 23	1
060.240.14*	1000	24 / 29	19 / 26	1
060.240.14A*	1000	24 / 29	24 / 29	1
060.240.15*	1000	29 / 32	14 / 23	1
060.240.16*	1000	29 / 32	19 / 26	1
060.240.17	1000	29 / 32	29 / 32	1
060.240.17A*	1000	34 / 35	24 / 29	1
060.240.18*	2000	24 / 29	19 / 26	1
060.240.18A*	2000	24 / 29	14 / 23	1
060.240.19*	2000	29 / 32	24 / 29	1
060.240.20	2000	29 / 32	29 / 32	1
060.240.20A*	2000	34 / 35	19 / 26	1
060.240.20B*	2000	34 / 35	24 / 29	1
060.240.26*	3000	34 / 35	19 / 26	1
060.240.26A*	3000	34 / 35	24 / 29	1
060.240.29*	5000	34 / 35	24 / 29	1
060.240.29B*	5000	34 / 35	19 / 26	1
060.240.32*	10000	34 / 35	24 / 29	1
060.240.33*	10000	45 / 40	24 / 29	1
060.240.33A*	10000	55 / 44	24 / 29	1
060.240.34*	20000	55 / 44	24 / 29	1



Buy Now!



PS : All 4 Neck Flasks can be provided in ASTM joint sizes on request.

These flasks can also be provided in other joint sizes

* Non DIN

Flasks

Flask, Round Bottom, Four Necks Parallel

- Complies with DIN 12392
- Centre Neck and Three Parallel Side Neck
- With Interchangeable Joints
- Made From Heat Resistant, Low Expansion 3.3 Borosilicate Glass



Item Code	Capacity (ml)	Interchangeable		Pack of
		Centre Neck	Side Neck	
061.240.00*	50	24 / 29	24 / 29	1
061.240.00B*	100	24 / 29	14 / 23	1
061.240.00BA*	100	24 / 29	19 / 26	1
061.240.01*	250	24 / 29	19 / 26	1
061.240.01A*	250	19 / 26	19 / 26	1
061.240.01B*	250	24 / 29	14 / 23	1
061.240.02*	250	29 / 32	14 / 23	1
061.240.02A*	250	29 / 32	29 / 32	1
061.240.03*	500	24 / 29	19 / 26	1
061.240.03A*	500	19 / 26	19 / 26	1
061.240.03B*	500	24 / 29	24 / 29	1
061.240.03C*	500	24 / 29	14 / 23	1
061.240.04*	500	29 / 32	14 / 23	1
061.240.05	500	29 / 32	19 / 26	1
061.240.05A	500	29 / 32	29 / 32	1
061.240.05B*	500	34 / 35	24 / 29	1
061.240.06*	1000	24 / 29	19 / 26	1
061.240.06A*	1000	24 / 29	14 / 23	1
061.240.06B*	1000	24 / 29	24 / 29	1
061.240.07*	1000	29 / 32	14 / 23	1
061.240.08*	1000	29 / 32	19 / 26	1
061.240.09*	1000	29 / 32	24 / 29	1
061.240.10	1000	29 / 32	29 / 32	1
061.240.10A*	1000	34 / 35	24 / 29	1
061.240.10B	1000	45 / 40	29 / 32	1
061.240.11*	2000	24 / 29	19 / 26	1
061.240.11A*	2000	24 / 29	14 / 23	1
061.240.12*	2000	29 / 32	19 / 26	1
061.240.13	2000	29 / 32	29 / 32	1
061.240.14*	2000	34 / 35	19 / 26	1
061.240.15*	2000	34 / 35	24 / 29	1
061.240.15B	2000	45 / 40	29 / 32	1
061.240.16B*	3000	29 / 32	29 / 32	1
061.240.17*	3000	34 / 35	19 / 26	1
061.240.18*	3000	34 / 35	24 / 29	1
061.240.18B*	3000	45 / 40	29 / 32	1
061.240.19A*	5000	34 / 35	19 / 26	1
061.240.19*	5000	34 / 35	24 / 29	1
061.240.26*	10000	34 / 35	24 / 29	1
061.240.26A*	10000	45 / 40	24 / 29	1
061.240.26B*	10000	55 / 44	24 / 29	1
061.240.27*	20000	55 / 44	24 / 29	1

PS : All 4 Neck Flasks can be provided in ASTM joint sizes on request.

* Non DIN

Диаэм, Москва ■ ул. Магаданская, д. 7, к. 3 ■ тел./факс: 8 (800) 234-0508 ■ sales@dia-m.ru

С.-Петербург
spb@dia-m.ru

Новосибирск
nsk@dia-m.ru

Воронеж
vrn@dia-m.ru

Йошкар-Ола
nba@dia-m.ru

Красноярск
krsk@dia-m.ru

Казань
kazan@dia-m.ru

Ростов-на-Дону
rnd@dia-m.ru

Екатеринбург
ekb@dia-m.ru

Кемерово
kemerovo@dia-m.ru

Нижний Новгород
nnovgorod@dia-m.ru

