

Vacuum Pump *New*

Salient Features:

- Valves made of SS 316 material Built-in micro suction filter
- All parts made from special graded Aluminum die cast material for light weight and good strength
- Totally oil-free construction
- No lubricant required
- Diaphragms are made of special 2-ply nylon reinforced neoprene rubber.
- Noiseless performance
- Extra large bearings for trouble-free and smooth running.
- Absolutely portable
- Practically maintenance –free



Item Code	Volt	Hz
970.DNEU.01	230	50



Buy Now!

Technical Specification :

Material	Aluminum Die Cast
Vacuum	22" Hg
Flow	15 (Ltrs/Min)
Pressure	25 Psig
Connection	10 mm Hose
Weight	2.5 Kg
Motor	1/20 hp
Dimension	150x110x140

000 «Диаэм»

Москва

ул. Магаданская, д. 7, к. 3 ■ тел./факс: (495) 745-0508 ■ sales@dia-m.ru

www.dia-m.ru

С.-Петербург
+7 (812) 372-6040
spb@dia-m.ru

Новосибирск
+7 (383) 328-0048
nsk@dia-m.ru

Воронеж
+7 (473) 232-4412
vrn@dia-m.ru

Йошкар-Ола
+7 (927) 880-3676
nba@dia-m.ru

Красноярск
+7 (923) 303-0152
krsk@dia-m.ru

Казань
+7 (843) 210-2080
kazan@dia-m.ru

Ростов-на-Дону
+7 (863) 303-5500
rnd@dia-m.ru

Екатеринбург
+7 (912) 658-7606
ekb@dia-m.ru

Кемерово
+7 (923) 158-6753
kemerovo@dia-m.ru

Армения
+7 (094) 01-0173
armenia@dia-m.ru

