

8-Channel Decapper

PRODUCT INTRODUCTION:

By using decapper, the opening of screw cap becomes more convenient. The head is arranged according to SBS standard to achieve that 8 caps in an individual line are opened / closed together. Torque limiting structure of decapper ensures caps sealing completely, but also avoids the situation of damaging sealing ring because of excessive torsion. The decapper is able to open /close Cryogenic Vials' caps in 3 seconds.

FEATURES OF PRODUCTS:

- * Open or close 8 vial's caps in 3 seconds.
- * Avoid touching samples or reagents. Reduce the risks of biochemical reagents and biological samples
- * It is the lightest among all handheld decappers. Alleviate fatigue caused by long-term use.
- * Fit for Biologix internal vials: 89-4233, 89-4253, 89-4283, 89-4213.

Finger rest / The ergonomically designed finger rest takes the load off the hand.

One hand operation Automatic tripping out the caps

CKS-4008



Internal rotation
external rotation caps available

