



Pharmacy Refrigerator

HLM-320

Scope of Application :

Suitable for the storage of biological products and items that need to be stored at 2° C ~8° C temperature. This product is applicable to pharmacies, epidemic prevention stations, and health centers as key examples.

Innovative Design

- ◉ Green Environmental Protection
- ◉ Digital Display with 0.1 Precision
- ◉ Large Capacity
- ◉ Small Foot Print
- ◉ Ergonomics Design
- ◉ Best in Class Performance

Product Advantages



Green Environmental Protection

The unit is designed and manufactured with a hydrocarbon compressor and hydrocarbon refrigerant resulting in an eco-friendlier refrigeration system



Large Capacity

323L large capacity storage space



Superior Performance

Indoor temperature is precisely controlled between 2°C to 8°C



Ergonomic Design




- Built-in lighting
- Door seals are removable
- Built-in door handle, no need for aftersales maintenance
- 4 adjustable shelves
- Professional and elegant appearance




Product Parts Diagram



Temperature Curve at 32°C Ambient Temperature



Product Dimensions



Technical Parameters

Model	Power Supply (V/Hz)	Power (W)	External Dimensions (W*D*H)(mm)	Internal Dimensions (W*D*H)(mm)	Effective Volume (L)	N.W/G.W (kg)	Shelf
HLM-320	220~240, 50	260	620*615*1838	532*445*1330	323	65/75	4

ДИАМ, Москва ■ ул. Магаданская, д. 7, к. 3 ■ тел./факс: 8 (800) 234-0508 ■ sales@dia-m.ru

ДИАМ
современная лаборатория

С.-Петербург
spb@dia-m.ru

Новосибирск
nsk@dia-m.ru

Воронеж
vrn@dia-m.ru

Йошкар-Ола
nba@dia-m.ru

Красноярск
krsk@dia-m.ru

Казань
kazan@dia-m.ru

Ростов-на-Дону
rnd@dia-m.ru

Екатеринбург
ekb@dia-m.ru

Кемерово
kemerovo@dia-m.ru

Нижний Новгород
nnovgorod@dia-m.ru

мобильное приложение



www.dia-m.ru